

# ORTHOPEDIC ATTACHMENT

MI 112



Hip Replacement



Tibia Nailing



Femur Nailing

## Objectives

Ortho Attachment allows the table to be used for the lower extremity procedures. Available at competitive prices our ortho attachment is purely made of stainless steel which resist corrosion.

## Key Specifications

- This Orthopedic Extension device can be used with any operating table as it independently stands on the floor.
- Specially Design for Intra Operative X-Ray photos & image intensifier set use.
- In addition, a pelvic support plate is supported on top of the frame and is formed with an asymmetrical construction. The pelvic support plate includes a cutout area on a front corner thereof permitting additional access to the affected leg of a patient. Removable Pelvic rest & X-Ray translucent pelvic stem.
- Stirring Instruments allows Flexion/Raising of Legs & Stirring/String ball joints locked from foot end.
- Two Traction Devices on sliding spot with ball joint allowing flexions & rotations of foot with traction shoe with metallic sole.
- Padded orthopedic knee Crutches/Provides support to the legs of the patient.
- Traction Unit with supporter is so easily & collapsible that it requires a very small place.
- The plate may be flipped over to position the cutout area on the opposite side of the patient, depending on the location of the affected leg.

---

**For Further Information, contact:**



487/67, 2<sup>nd</sup> Floor, Glass Street,  
National Market, Peeragarhi,  
New Delhi, India – 110087

Tele: 91-11-25255415 / 416

Fax: 91-11-25255417

Mail: [info@meditech-india.com](mailto:info@meditech-india.com)

Web: [www.meditechindia.in](http://www.meditechindia.in)

**Online assistance:**



MSN: marketing\_meditech



Skype: exportmeditech