



## TECHNICAL DATA

	Z 500/Z 500 M
Light Intensity @ 1 m	150000 Lux
Number of LED's	90
Life-Span of LED's	50000 hrs.
Focus Dia	15-25 cm
Light Head Size	62.5 cm
Depth of Illumination	120 cm
Color Rendering Index Ra (1)	92
Color Temperature	4300K
Illuminating Adjustment	10% - 100%
Height Adjustment	120 cm
Total Power Consumption	50 W
Focusing	Adjustable
Power Supply	160 V - 240 V DC, 50 Hz
Control Unit	SMPS